



ECG

SUCTION CHEST ELECTRODE



UMDNS

11439

Description

The suction chest electrodes are characterized by the maximum hold on the skin of the patient and the best performance in terms of registration and duration. The connector is universal and it is available in 3 different sizes.

Different configurations available on request

Materials

Sensor AgAgCl, Ni

Ball pump PVC

Packaging LDPE

Safety

Biocompatibility
(UNI EN ISO 10993) Yes

Latex No

Phthalates No

Specifications

Generic code
SUXXXXXX,
BALLXXXXXX

Electrode diameter
AgAgCl from 15 to 30 mm,
Ni 24 mm

Electrode colour ●

Connector Universal connection

Environment

RoHS II ✓

REACH ✓