



EEG

CASTED AUTOCLAVABLE CUP ELECTRODE



UMDNS

11440

Description

The Casted version of electrodes cup is characterized by the durability of the product over time.

The electrode is waterproof, its special shape protects the electrode-lead wire connection from possible contacts with water. It is also autoclavable, provided with a silicone ultraflexible cable.

- Its special shape protects the electrode-lead wire connection from possible contacts with water ensuring a longer duration of the product;
- Autoclavable electrode, equipped with silicone ultraflexible lead wire;
- Pack of 10 colours.


Materials

Electrode	AgAgCl - Ag - Au
Lead wire	Cu/Sn with SI jacket
Connector	CuZn/Au with PP jacket
Packaging	PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	DXXXXXXXXXX
Electrode diameter	10 mm
Cable length	From 100 to 250 cm
Cable colour	
Connector	Touch proof DIN 42802

Environment

RoHS II	✓
REACH	✓

Different configurations available on request