



EMG

GROUND BAND ELECTRODE



UMDNS

11441

Description

The ground electrode has an effective role in order to have the best quality recording.

Spes Medica offers two different kinds of ground bracelet: the version for the patient's best comfort and the pure steel version that can be used dry.

In every package you can find 1 meter connection lead wire which connects the band to the headbox.

Materials

Sensor	Stainless Steel
Lead wire	Cu/Sn with PVC jacket
Connector	PE
Packaging	PP

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Latex	No
Nickel	No
Phthalates	No

Specifications

Generic code	GDRXXXXXX
Size	From 28 to 60 cm
Connector	Touch proof DIN 42802

Environment

RoHS II	✓
REACH	✓

Other configurations available on request