



EEG



PSG

## SLEEPCAP

### Prewired electrode cap



UMDNS

17554

### Description

Thanks to Sleepcap is possible to have the advantages of a pre-wired cap even in sleep recording.

The reason is that the electrode and the cap are totally made by soft material for the best patient's comfort and for sleep recordings.

- Net and electrodes completely made by soft material for the best comfort for the patient;
- Configuration with inside or outside wires;
- On request available to 256 electrode even in the smaller sizes;
- Available with different type of connectors and configuration for the most EEG equipments.

### Materials

Electrode	Tin (on request AgAgCl, sintered AgAgCl or Gold)
Electrode holder	PU
Lead wire	Cu/Sn with PVC jacket
Headcap net	PU
Box	PP

### Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

### Specifications

Generic code	CSSXXXXXXXXX
Electrode number	From 9 to 256
Connector	Available for every type of connectors

### Environment

RoHS II	✓
REACH	✓

Different configurations available on request

