



IOM

DRAGONFLY LARYNGEAL ELECTRODE



UMDNS

11441

Description

Spes Medica, thanks to the collaboration with Neurovision Medical presents a complete range of laryngeal electrodes. The electrodes can be applied in the whole range of endotracheal tubes, the electrode material is soft to avoid the risk to damage the vocal cords.

The electrodes are integrated with a long disposable twisted lead wire with Touch proof connector to fit in the most used EMG units.

- Silver sensor;
- Positioning marker.

Materials

Electrode	Ag
Substrate	PET
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PVC jacket
Packaging	Tyvek/PET-PE

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Sterilization method	Ehtylene oxide (EtO)
Latex	No
Phthalates	No

Specifications

Generic codes	LSEXXXXXXXXXXX
Number of channels	1 - 2
Height	From 34 to 42 mm
Width	From 46 to 63 mm
Cable length	75 - 200 cm
Connector	Touch proof DIN 42802

Cable colour ● ●

Environment

RoHS II	✓
REACH	✓

Different configurations available on request