



EEG



IOM

DISPOSABLE SUBDERMAL BENT NEEDLE ELECTRODES



UMDNS

11440

Description

The needle hook is especially used in situations where you need to ensure the positioning of the same for long recordings such as IOM or brain death.

The product has been recently improved by automating the production and redesigning the bending angle, ensuring better contact with the skin and lower impedance of the electrode contact.

Materials

Needle	Stainless steel
Needle holder	PVDF
Needle cover	PVC
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Au with PVC jacket
Packaging	PP

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

Specifications

Generic code	MNXXXXXXXXXX
Needle diameter	From 0,20 to 0,40 mm
Needle length	From 10 to 22 mm
Cable length	From 8 to 250 cm
Cable colour	● ● ● ● ● ●
Connector	Touch proof DIN 42802, 1 mm Pin

Environment

RoHS II	✓
REACH	✓

Different configurations available on request