



EEG



IOM

NEURODART: DISPOSABLE SUBDERMAL NEEDLE ELECTRODE



UMDNS

11440

Description

Spes Medica offers the quality you expect and help deliver the patient comfort that you need.

The needle electrode material is biocompatible stainless steel: the best compromise between clear recording with low noise and cost effective.

With a variety of needle lengths and types, available in several lead wire lengths, we have the right set of needle electrodes for your application needs.

- Made of high-tensile strength stainless steel;
- Cleared for both stimulating and recording;
- High clear recording due to the low impedance of contact;
- Best sharpness of the needle for the best placement.

Materials

Needle	Stainless Steel
Needle holder	PVDF
Needle cover	PVC
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn/Au with PVC jacket
Packaging	Medical paper/PET

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxyde (EtO)
Latex	No
Phthalates	No

Specifications

Generic code	MNXXXXXXXXXX
Needle diameter	From 0,20 to 0,40 mm
Needle length	From 10 to 22 mm
Cable length	From 8 to 250 cm
Cable colour	
Connector	Touch proof DIN 42802

Environment

RoHS II	✓
REACH	✓

Different configurations available on request