



EMG

GROUND FLAT ELECTRODE



UMDNS

11441

Description

These 30 mm diameter electrodes ensure the maximum flexibility to the operator since they can be positioned all over the patient's body.

Materials

Sensor	Ni
Needle handgrip	PE
Lead wire	Cu/Sn with PVC jacket
Connector	PE
Packaging	PP

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Latex	No
Nickel	No
Phthalates	No

Specifications

Generic code	GDRXXXXXXXX
Electrode diameter	30 mm
Lead wire length	125 cm
Connector	Touch proof DIN 42802

Environment

RoHS II	✓
REACH	✓

Other configurations available on request