



EEG

MEDCAP

Prewired electrode cap



UMDNS

17554

Description

The main features of this cap is the dependability tested in many years with all different type of EEG equipment, the possibility to repair the cap quickly and easily and the great flexibility of this type of product for the EEG lab.

- Elastic material headcap;
- Easy colour coding system for different sizes;
- Easy to wash and replace the electrode.

Materials

Electrode	Tin (on request AgAgCl, sintered AgAgCl or Gold)
Electrode holder	LDPE
Lead wire	PVC
Headcap net	PU
Packaging	PP

Safety

Biocompatibility <i>(UNI EN ISO 10993)</i>	Yes
Latex	No
Nickel	No
Phthalates	No

Specifications


Generic code	CAMXXXXXXXX
Electrode number	From 9 to 256
Connector	Available for every type of connector

Environment


RoHS II	✓
REACH	✓

Different configurations available on request

Adult



Infant



Neonatal

