



EEG



IOM

## TWISTED NEEDLE ELECTRODES



UMDNS

11441

### Description

Disposable twisted monopolar subdermal needles.  
Available in set of 2, 5 and 10 twisted needles.  
The twisted needle electrodes are easy to recognize thank to the colour coding of the lead wire.

- Special sharpening for a painless and atraumatic needle insertion;
- 100% quality control of the product;
- Very thin shrinking tube as alternative to the needles with moulded hub;
- Pair of two needles twisted with cables and different combinations of colours.

*Different configurations available on request*

### Materials

Needle	Stainless Steel
Needle handgrip	PVDF
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PVC jacket
Packaging	Medical paper/PET-PE

### Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

### Specifications

Generic codes	MNXXXXXXXXXX, TTXXMNXXXXD
Needle diameter	From 0,20 to 0,40 mm
Needle lenght	From 10 to 22 mm
Cable length	From 100 to 250 cm
Cable colour	
Connector	Touch Proof DIN 42802

### Environment

RoHS II	✓
REACH	✓