




EMG

PUDENDAL ELECTRODE



UMDNS

11441

Description

The electrode for the pelvic floor can be used for both stimulation and recording. It is an invasive device, previously placed on the examination glove, used by the operator on patients with pelvic floor diseases.

- Possibility to stimulate and record the pelvic floor;
- Biocompatible materials, latex free.

Materials

Electrode	Ag
Support	PET
Packaging	Medical paper/PP

Safety

Biocompatibility (UNI EN ISO 10993)	Yes
Latex	No
Phthalates	No

Specifications

Generic code	PUDXXXX
Cable length	2 m
Connector	Touch proof DIN 42802, 5 pole DIN 45322
Connection cable	2001PUDXXXXXX

Environment

PVC-Free	✓
RoHS II	✓
REACH	✓

Different configurations available on request